

OpenSesame

A graphical experiment builder for the social sciences

Binus University Workshop 2022

Sebastiaan Mathôt

https://osdoc.cogsci.nl/binus2022

Sponsored by





A round of introduction

About me

- Sebastiaan Mathôt
- Originally from Amsterdam (Netherlands)
- Lived for a while in Marseille (France)
- Now working in Groningen (Netherlands)

About you

Menti.com → 20 94 008



Workshop schedule



Workshop schedule

Four days

13:00 - 17:00

April 1 (Today): Introduction

April 20: Online experiments

April 21: Eye tracking

April 27: Data Analysis

All days: working on your own experiment!



Workshop schedule

https://osdoc.cogsci.nl/binus2022

Binus University 2022 workshop

Edit on GitHub

- Practical information
- Description
- . Day 1: Introduction (April 1)
- Day 2: Online experiments (April 20)
- Day 3: Eye-tracking experiments (April 21)
- Day 4: Data analysis (April 27)
- Suggested experiments

Practical information

- · Host: Binus University
- Location: online
- · Dates:
 - o Day 1: Friday, April 1st, 2022, 08:00-12:00 CEST, 13:00 17:00 Jakarta
 - o Day 2: Wednesday, April 20th, 2022, 08:00-12:00 CEST, 13:00 17:00 Jakarta
 - o Day 3: Thursday, April 21st, 2022, 08:00-12:00 CEST, 13:00 17:00 Jakarta
 - Day 4: Wednesday, April 27th, 2022, 08:00-12:00 CEST, 13:00 17:00 Jakarta
- · Presenter: Sebastiaan Mathôt
- · Spreadsheet with participant overview

Description

In this four-day, hands-on, online workshop, you will learn how to implement psychological experiments with the open-source software OpenSesame. You will learn:

- · How to run experiments online as well as in a traditional laboratory set-up.
- · The limitations and advantages of online and laboratory-based experiments.
- How to include eye tracking in laboratory-based experiments. (And a sneak-peak at eye-tracking in online experiments!)
- · How to analyze data collected from online and laboratory-based experiments.

Finally, using the skills that you will learn during the workshop, you will design and implement an experiment for your own research, of course with assistance from us! For this purpose, please already think about what kind of experiment you'd like to create, and






A graphical experiment builder

Drag-and-drop, point-and-click

Complement with scripting

- Python for lab experiments
- JavaScript for online experiments





Open source


Free of charge





Cross-platform

- Windows
- Mac OS
- Linux
- Browser (runtime only)





Fits many types of research

- Psychophysics
- Neurophysiology
- Neuroimaging
- Social psychology
- Clinical applications
- Eye tracking
- Pupillometry





Support




Support

Documentation

https://osdoc.cogsci.nl/

Community forum

- https://forum.cogsci.nl/
- Shared with JASP, PyGaze, JATOS, and other free software
- 9,000+ members, daily activity





Using OpenSesame



Items

An experiment is a collection of items


Items are building blocks

One item does one thing

- sketchpad → visual stimuli
- keyboard_response → records key presses
- logger → logs data
- loop → repeats another item while varying independent variables
- sequence → runs a number of other items sequentially
- etc.



Items






Items

Items have a graphical interface and a script

- A simple language ("OpenSesame script")
- Not Python!

Switch between graphical interface and script





Variables can be defined

- In a loop item (independent variables)
- In a keyboard_response (response time, accuracy, etc.)

You can refer to these later

- For example in a sketchpad item
- Using [variable_name] notation


This adds considerable flexibility





	gaze_cue	target_letter	target_pos	congruency	correct_response	dist_pos
1	left	F	-300	congruent	Z	300
2	right	F	-300	incongruent	Z	300
3	left	Н	-300	congruent	m	300
4	right	Н	-300	incongruent	m	300
5	left	F	300	incongruent	Z	-300
6	right	F	300	congruent	Z	-300
7	left	Н	300	incongruent	m	-300
8	right	Н	300	congruent	m	-300









▼ Flush pending key presses at sequence start

Item name	Run if	
■ trial_sequence		
fixation_dot	always	
neutral_gaze	always	
gaze_cue	always	
target	always	
keyboard_response	always	
incorrect_sound	[correct] = 0	
🐷 logger	always	





Let's get to work!