



Eye tracking Binus University Workshop Day 3

https://osdoc.cogsci.nl/binus2022


Sebastiaan Mathôt











Today (day 3)




- Before the break
 - A general introduction to eye tracking
 - Working on a self-paced reading task
- After the break
 - Gaze-contingent self-paced reading
 - Continue working on your own experiment

Eye tracking




- What do eye trackers record?
- Always: gaze position (time series)
 - x, y coordinates
 - Relative to a display
 - Or relative to a world camera
 - Or uncalibrated
- Usually: pupil size (time series)
 - In arbitrary units
 - Or in millimeters of diameter
- Usually: an event log
 - Timestamped messages



Eye tracking



- Eye tracking is not standardized
 - Different log-file formats
 - Different sampling rates
 - Different software-development kits (SDKs)
 - Technical idiosyncrasies
 - No standard analysis pipeline
- Attempts at standardization
 - Building experiments
 - PyGaze/ OpenSesame
 - iohub/ PsychoPy
 - Analyzing data (script-based)
 - gazeR
 - python-eyelinkparser



Eye tracking



- An eye-tracker setup consists of
 - Experiment PC: runs the experiment
 - Recording PC: connected to eye tracker
 - Bi-directional communication
- For some eye trackers the experiment and recording computers are the same




EyeLink 1000 setup Source: SR Research

PyGaze



- PyGaze is a Python library for eye tracking
 - Supports different brands
 - Integrates with OpenSesame
- We will use PyGaze to implement a fancy gazecontingent, self-paced reading experiment
 - But we'll start slow



Eye tracking in a browser




- High-quality eye tracking requires a laboratory setup
- But poor quality eye tracking can be done in a browser
 - Using the webcam
 - WebGazer.js is the bestknown implementation
 - Works with OSWeb [1]



Hands-on workshop



- As said: We'll implement a self-paced reading task
- I'll suggest some pointers
- But of course the details are up to you!




Source: MagicFlyingBunnies at deviantart

Hands-on workshop



Independent variable: sentence complexity (?)




Dependent variable: response time

Hands-on workshop



- Start from the Default template
- At the start of the experiment
 - pygaze_init → initialize the eye tracker in advanced dummy mode
- At the start of each trial
 - pygaze_drift_correct
 - pygaze_start_recording
- At the end of each trial
 - pygaze_log
 - pygaze_stop_recording



Source: MagicFlyingBunnies at deviantart



Let's get to work!

Slides: https://osdoc.cogsci.nl/binus2022

Adding gaze contingency



