



university of
 groningen

faculty of behavioural
and social sciences



Eye tracking

Leuphana workshop Day 3

<https://osdoc.cogsci.nl/leuphana2021>

Sebastiaan Mathôt



cogsci.nl/smathot



[cognitivescience](https://www.facebook.com/cognitivescience)



[@smathot](https://github.com/smathot)



[@cogscinl](https://twitter.com/cogscinl)



[sebastiaanmathot](https://www.youtube.com/sebastiaanmathot)

Today (day 3)




- Before the break
 - A general introduction to eye tracking
 - Working on a self-paced reading task
- After the break
 - Gaze-contingent self-paced reading
 - Continue working on your own experiment
 - Q&A
 - Workshop end!

Eye tracking



- What do eye trackers record?
- Always: gaze position (time series)
 - x, y coordinates
 - Relative to a display
 - Or relative to a world camera
 - Or uncalibrated
- Usually: pupil size (time series)
 - In arbitrary units
 - Or in millimeters of diameter
- Usually: an event log
 - Timestamped messages



Eye tracking

- Eye tracking is not standardized
 - Different log-file formats
 - Different sampling rates
 - Different software-development kits (SDKs)
 - Technical idiosyncrasies
 - No standard analysis pipeline
- Attempts at standardization
 - Building experiments
 - PyGaze/ OpenSesame
 - iohub/ PsychoPy
 - Analyzing data (script-based)
 - gazeR
 - python-eyelinkparser

Eye tracking



- An eye-tracker setup consists of
 - Experiment PC: runs the experiment
 - Recording PC: connected to eye tracker
 - Bi-directional communication
- For some eye trackers the experiment and recording computers are the same




EyeLink 1000 setup
Source: SR Research

PyGaze




- PyGaze is a Python library for eye tracking
 - Supports different brands
 - Integrates with OpenSesame
- We will use PyGaze to implement a fancy gaze-contingent, self-paced reading experiment
 - But we'll start slow



Hands-on workshop



- As said: We'll implement a self-paced reading task
- I'll suggest some pointers
- But of course the details are up to you!




Source: MagicFlyingBunnies at deviantart

Hands-on workshop



Independent variable: sentence complexity (?)




Dependent variable: response time

Hands-on workshop



- Start from the *Default template*
- At the start of the experiment
 - **pygaze_init** → initialize the eye tracker in advanced dummy mode
- At the start of each trial
 - **pygaze_drift_correct**
 - **pygaze_start_recording**
- At the end of each trial
 - **pygaze_log**
 - **pygaze_stop_recording**



Source: MagicFlyingBunnies at deviantart



Let's get to work!

Slides: <https://osdoc.cogsci.nl/leuphana2021>