



university of  
 groningen

faculty of behavioural  
 and social sciences



# Eye tracking

## Leuphana workshop Day 3

<https://osdoc.cogsci.nl/leuphana2021>

Sebastiaan Mathôt



[cogsci.nl/smathot](https://cogsci.nl/smathot)



[cognitivescience](https://www.facebook.com/cognitivescience)



[@smathot](https://github.com/smathot)



[@cogscinl](https://twitter.com/cogscinl)



[sebastiaanmathot](https://www.youtube.com/sebastiaanmathot)

# Today (day 3)




- Before the break
  - A general introduction to eye tracking
  - Working on a self-paced reading task
- After the break
  - Gaze-contingent self-paced reading
  - Continue working on your own experiment
  - Q&A
  - Workshop end!

# Eye tracking




- What do eye trackers record?
- Always: gaze position (time series)
  - x, y coordinates
  - Relative to a display
  - Or relative to a world camera
  - Or uncalibrated
- Usually: pupil size (time series)
  - In arbitrary units
  - Or in millimeters of diameter
- Usually: an event log
  - Timestamped messages



# Eye tracking



- Eye tracking is not standardized
  - Different log-file formats
  - Different sampling rates
  - Different software-development kits (SDKs)
  - Technical idiosyncrasies
  - No standard analysis pipeline
- Attempts at standardization
  - Building experiments
    - PyGaze/ OpenSesame
    - iohub/ PsychoPy
  - Analyzing data (script-based)
    - gazeR
    - python-eyelinkparser



# Eye tracking



- An eye-tracker setup consists of
  - Experiment PC: runs the experiment
  - Recording PC: connected to eye tracker
  - Bi-directional communication
- For some eye trackers the experiment and recording computers are the same




EyeLink 1000 setup  
Source: SR Research

# PyGaze



- PyGaze is a Python library for eye tracking
  - Supports different brands
  - Integrates with OpenSesame
- We will use PyGaze to implement a fancy gaze-contingent, self-paced reading experiment
  - But we'll start slow






# Hands-on workshop



- As said: We'll implement a self-paced reading task
- I'll suggest some pointers
- But of course the details are up to you!




Source: MagicFlyingBunnies at deviantart

# Hands-on workshop



Independent variable: sentence complexity (?)




Dependent variable: response time



# Hands-on workshop



- Start from the *Default template*
- At the start of the experiment
  - **pygaze\_init** → initialize the eye tracker in advanced dummy mode
- At the start of each trial
  - **pygaze\_drift\_correct**
  - **pygaze\_start\_recording**
- At the end of each trial
  - **pygaze\_log**
  - **pygaze\_stop\_recording**



Source: MagicFlyingBunnies at deviantart



# Let's get to work!

Slides: <https://osdoc.cogsci.nl/leuphana2021>

# Exercise

